

Recombinant Human DCTN6

Catalog No: #GP11449



Package Size: #GP11449-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DCTN6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	dynactin 6
Other Names	WS3; p27; WS-3
Accession No.	Swissprot:O00399Gene Accession:BC013175
Uniprot	O00399
GeneID	10671;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene contains an RGD (Arg-Gly-Asp) motif in the N-terminal region, which confers adhesive properties to macromolecular proteins like fibronectin. It shares a high degree of sequence similarity with the mouse homolog, which has been suggested to play a role in mitochondrial biogenesis. The exact biological function of this gene is not known.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only