

Recombinant Human DDX53

Catalog No: #GP11454



Package Size: #GP11454-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DDX53
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 253-428 amino acids of human DEAD (Asp-Glu-Ala-Asp) box polypeptide 53
Target Name	DEAD (Asp-Glu-Ala-Asp) box polypeptide 53
Other Names	CAGE; CT26
Accession No.	Swissprot:Q86TM3Gene Accession:BC067878
Uniprot	Q86TM3
GeneID	168400;
Storage	-20~-80°C, pH 7.6 PBS

Background

This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only