

Recombinant Human DHRS3

Catalog No: #GP11462



Package Size: #GP11462-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DHRS3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 30-169 amino acids of human dehydrogenase/reductase (SDR family) member 3
Target Name	dehydrogenase/reductase (SDR family) member 3
Other Names	SDR1; RDH17; Rsd1; DD83.1; SDR16C1; retSDR1
Accession No.	Swissprot:O75911Gene Accession:BC002730
Uniprot	O75911
GeneID	9249;
Storage	-20~-80°C, pH 7.6 PBS

Background

Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only