

Recombinant Human DNAJB4

Catalog No: #GP11469



Package Size: #GP11469-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DNAJB4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	DnaJ (Hsp40) homolog, subfamily B, member 4
Other Names	DjB4; HLJ1; DNAJW
Accession No.	Swissprot:Q9UDY4Gene Accession:BC034721
Uniprot	Q9UDY4
GeneID	11080;
Storage	-20~-80°C, pH 7.6 PBS

Background

DNAJB4 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only