

Recombinant Human EDIL3

Catalog No: #GP11489



Package Size: #GP11489-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EDIL3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human EGF-like repeats and discoidin I-like domains 3
Target Name	EGF-like repeats and discoidin I-like domains 3
Other Names	DEL1
Accession No.	Swissprot:O43854Gene Accession:BC030828
Uniprot	O43854
GeneID	10085;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only