

Recombinant Human EIF5

Catalog No: #GP11496



Package Size: #GP11496-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EIF5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 233-431 amino acids of human eukaryotic translation initiation factor 5
Target Name	eukaryotic translation initiation factor 5
Other Names	EIF-5; EIF-5A
Accession No.	Swissprot:P55010Gene Accession:BC007728
Uniprot	P55010
GenID	1983;
Storage	-20~-80°C, pH 7.6 PBS

Background

Eukaryotic translation initiation factor-5 (EIF5) interacts with the 40S initiation complex to promote hydrolysis of bound GTP with concomitant joining of the 60S ribosomal subunit to the 40S initiation complex. The resulting functional 80S ribosomal initiation complex is then active in peptidyl transfer and chain elongations

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only