

Recombinant Human KIAA1324

Catalog No: #GP11498



Package Size: #GP11498-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KIAA1324
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 425-674 amino acids of human KIAA1324
Target Name	KIAA1324
Other Names	EIG121
Accession No.	Swissprot:Q6UXG2Gene Accession:NP_065826
Uniprot	Q6UXG2
GeneID	57535;
Storage	-20~-80°C, pH 7.6 PBS

Background

Expression of this gene is induced by estrogen and the encoded protein has been characterized as a transmembrane protein. The encoded protein has been found in to correlate with survival in certain carcinomas (PMID: 21102415) and may be important for cellular response to stress (PMID: 21072319). Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only