

Recombinant Human ELF1

Catalog No: #GP11502



Package Size: #GP11502-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ELF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 370-619 amino acids of human ELF1
Target Name	E74-like factor 1 (ets domain transcription factor)
Accession No.	Swissprot:P32519Gene Accession:BC030507
Uniprot	P32519
GeneID	1997;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only