

# Recombinant Human EDN1

Catalog No: #GP11517



Package Size: #GP11517-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human EDN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 74-212 amino acids of human endothelin 1
Target Name	endothelin 1
Other Names	ET1; PPET1; HDLCQ7
Accession No.	Swissprot:P05305Gene Accession:BC009720
Uniprot	P05305
GeneID	1906;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is proteolytically processed to release a secreted peptide termed endothelin 1. This peptide is a potent vasoconstrictor and is produced by vascular endothelial cells. Endothelin 1 also can affect the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only