

UBE2Z Antibody

Catalog No: #42999

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBE2Z Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total UBE2Z protein.
Immunogen Description	Fusion protein of human UBE2Z
Target Name	UBE2Z
Other Names	USE1; HOYS7
Accession No.	Swiss-Prot#: Q9H832 Gene ID: 65264
Uniprot	Q9H832
GeneID	65264;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #42999 at dilution 1/100.



Immunohistochemical analysis of paraffin-embedded Human Liver cancer tissue using #42999 at dilution 1/100.

Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

Note: This product is for in vitro research use only