

Recombinant Human ETV7

Catalog No: #GP11523



Package Size: #GP11523-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ETV7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-151 amino acids of human ETV7
Target Name	ets variant 7
Other Names	TEL2; TELB; TEL-2
Accession No.	Swissprot:Q9Y603Gene Accession:BC035853
Uniprot	Q9Y603
GeneID	51513;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the ETS family of transcription factors, which is a large group of evolutionarily conserved transcriptional regulators that play an important role in a variety of cellular processes throughout development and differentiation, and are involved in oncogenesis as well. This protein is predominantly expressed in hematopoietic tissues. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only