

Recombinant Human EZR

Catalog No: #GP11525



Package Size: #GP11525-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EZR
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-295 amino acids of human ezrin
Target Name	ezrin
Other Names	CVL; CVIL; VIL2; HEL-S-105
Accession No.	Swissprot:P15311Gene Accession:BC013903
Uniprot	P15311
GeneID	7430;
Storage	-20~-80°C, pH 7.6 PBS

Background

The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only