

# Recombinant Human FAAH2

Catalog No: #GP11526



Package Size: #GP11526-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FAAH2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human fatty acid amide hydrolase 2
Target Name	fatty acid amide hydrolase 2
Other Names	AMDD
Accession No.	Swissprot:Q6GMR7 Gene Accession:BC073922
Uniprot	Q6GMR7
GeneID	158584;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a fatty acid amide hydrolase that shares a conserved protein motif with the amidase signature family of enzymes. The encoded enzyme is able to catalyze the hydrolysis of a broad range of bioactive lipids, including those from the three main classes of fatty acid amides; N-acylethanolamines, fatty acid primary amides and N-acyl amino acids. This enzyme has a preference for monounsaturated acyl chains as a substrate.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only