

Recombinant Human FABP6

Catalog No: #GP11528



Package Size: #GP11528-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FABP6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	fatty acid binding protein 6, ileal
Other Names	ILBP; I-15P; I-BAP; ILBP3; ILLBP; I-BABP; I-BALB
Accession No.	Swissprot:P51161Gene Accession:BC022489
Uniprot	P51161
GeneID	2172;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the ileal fatty acid binding protein. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABP6 and FABP1 (the liver fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism. Transcript variants generated by alternate transcription promoters and/or alternate splicing have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only