

Recombinant Human FAP

Catalog No: #GP11538



Package Size: #GP11538-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FAP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human fibroblast activation protein, alpha
Target Name	fibroblast activation protein, alpha
Other Names	FAPA; DPPIV
Accession No.	Swissprot:Q12884Gene Accession:BC026250
Uniprot	Q12884
GeneID	2191;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a homodimeric integral membrane gelatinase belonging to the serine protease family. It is selectively expressed in reactive stromal fibroblasts of epithelial cancers, granulation tissue of healing wounds, and malignant cells of bone and soft tissue sarcomas. This protein is thought to be involved in the control of fibroblast growth or epithelial-mesenchymal interactions during development, tissue repair, and epithelial carcinogenesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only