

Recombinant Human FASTKD2

Catalog No: #GP11543



Package Size: #GP11543-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FASTKD2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 300-580 amino acids of human FAST kinase domains 2
Target Name	FAST kinase domains 2
Other Names	KIAA0971
Accession No.	Swissprot:Q9NYY8Gene Accession:BC001544
Uniprot	Q9NYY8
GeneID	22868;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only