

Recombinant Human FGL2

Catalog No: #GP11557



Package Size: #GP11557-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FGL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 204-436 amino acids of human FGL2
Target Name	fibrinogen like 2
Other Names	T49; pT49
Accession No.	Swissprot:Q14314Gene Accession:BC033820
Uniprot	Q14314
GeneID	10875;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a secreted protein that is similar to the beta- and gamma-chains of fibrinogen. The carboxyl-terminus of the encoded protein consists of the fibrinogen-related domains (FRED). The encoded protein forms a tetrameric complex which is stabilized by interchain disulfide bonds. This protein may play a role in physiologic functions at mucosal sites.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only