

Recombinant Human FHL2

Catalog No: #GP11560



Package Size: #GP11560-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FHL2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	four and a half LIM domains 2
Other Names	DRAL; AAG11; FHL-2; SLIM3; SLIM-3
Accession No.	Swissprot:Q14192Gene Accession:NP_963849
Uniprot	Q14192
GeneID	2274;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the four-and-a-half-LIM-only protein family. Family members contain two highly conserved, tandemly arranged, zinc finger domains with four highly conserved cysteines binding a zinc atom in each zinc finger. This protein is thought to have a role in the assembly of extracellular membranes. Also, this gene is down-regulated during transformation of normal myoblasts to rhabdomyosarcoma cells and the encoded protein may function as a link between presenilin-2 and an intracellular signaling pathway. Multiple alternatively spliced variants, encoding the same protein, have been identified.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only