

Recombinant Human FIS1

Catalog No: #GP11564



Package Size: #GP11564-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FIS1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 122 amino acids of human fission 1 (mitochondrial outer membrane) homolog (<i>S. cerevisiae</i>)
Target Name	fission 1 (mitochondrial outer membrane) homolog (<i>S. cerevisiae</i>)
Other Names	TTC11; CGI-135
Accession No.	Swissprot:Q9Y3D6Gene Accession:BC003540
Uniprot	Q9Y3D6
GenID	51024;
Storage	-20~-80°C, pH 7.6 PBS

Background

The balance between fission and fusion regulates the morphology of mitochondria. TTC11 is a component of a mitochondrial complex that promotes mitochondrial fission.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only