

Recombinant Human FMNL1

Catalog No: #GP11568



Package Size: #GP11568-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FMNL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human formin-like 1
Target Name	formin-like 1
Other Names	FMNL; FHOD4; KW-13; C17orf1; C17orf1B
Accession No.	Swissprot:O95466Gene Accession:BC021906
Uniprot	O95466
GeneID	752;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. An alternative splice variant has been described but its full length sequence has not been determined.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only