

Recombinant Human FRK

Catalog No: #GP11570



Package Size: #GP11570-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FRK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 42-208 amino acids of human fyn-related kinase
Target Name	fyn-related kinase
Other Names	GTK; RAK; PTK5
Accession No.	Swissprot:P42685Gene Accession:BC012916
Uniprot	P42685
GeneID	2444;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only