

Recombinant Human FSTL1

Catalog No: #GP11571



Package Size: #GP11571-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human FSTL1 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 86-308 amino acids of human follistatin-like 1 |
| Target Name | follistatin-like 1 |
| Other Names | FRP; FSL1; MIR198 |
| Accession No. | Swissprot:Q12841Gene Accession:BC000055 |
| Uniprot | Q12841 |
| GeneID | 11167; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only