

Recombinant Human GALK1

Catalog No: #GP11580



Package Size: #GP11580-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GALK1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human galactokinase 1
Target Name	galactokinase 1
Other Names	GK1; GALK; HEL-S-19
Accession No.	Swissprot:P51570Gene Accession:BC001166
Uniprot	P51570
GeneID	2584;
Storage	-20~-80°C, pH 7.6 PBS

Background

Galactokinase is a major enzyme for the metabolism of galactose and its deficiency causes congenital cataracts during infancy and presenile cataracts in the adult population.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only