

Recombinant Human PAEP

Catalog No: #GP11602



Package Size: #GP11602-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAEP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-180 amino acids of human progesterone-associated endometrial protein
Target Name	progesterone-associated endometrial protein
Other Names	GD; GdA; GdF; GdS; PEP; PAEG; PP14
Accession No.	Swissprot:P09466Gene Accession:BC069451
Uniprot	P09466
GenID	5047;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the lipocalin superfamily whose members share relatively low sequence similarity but have highly conserved exon/intron structure and three-dimensional protein folding. Most lipocalins are clustered on the long arm of chromosome 9. The encoded glycoprotein has been previously referred to as pregnancy-associated endometrial alpha-2-globulin, placental protein 14, and glycodelin, but has been officially named progesterone-associated endometrial protein. Three distinct forms, with identical protein backbones but different glycosylation profiles, are found in amniotic fluid, follicular fluid and seminal plasma of the reproductive system.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only