

Recombinant Human GNG2

Catalog No: #GP11607



Package Size: #GP11607-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GNG2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	guanine nucleotide binding protein (G protein), gamma 2
Accession No.	Swissprot:P59768Gene Accession:BC020774
Uniprot	P59768
GeneID	54331;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes one of the gamma subunits of a guanine nucleotide-binding protein. Such proteins are involved in signaling mechanisms across membranes. Various subunits forms heterodimers which then interact with the different signal molecules.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only