

# Recombinant Human MB

Catalog No: #GP11611



Package Size: #GP11611-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MB
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	myoglobin
Other Names	PVALB; myoglobin
Accession No.	Swissprot:P02144Gene Accession:BC014547
Uniprot	P02144
GeneID	4151;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the globin superfamily and is expressed in skeletal and cardiac muscles. The encoded protein is a haemoprotein contributing to intracellular oxygen storage and transcellular facilitated diffusion of oxygen. At least three alternatively spliced transcript variants encoding the same protein have been reported.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only