

Recombinant Human LAMTOR5

Catalog No: #GP11618



Package Size: #GP11618-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LAMTOR5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	late endosomal/lysosomal adaptor, MAPK and MTOR activator 5
Other Names	XIP; HBXIP
Accession No.	Swissprot:O43504Gene Accession:BC062619
Uniprot	O43504
GeneID	10542;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that specifically complexes with the C-terminus of hepatitis B virus X protein (HBx). The function of this protein is to negatively regulate HBx activity and thus to alter the replication life cycle of the virus.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only