

Recombinant Human HOMER1

Catalog No: #GP11630



Package Size: #GP11630-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HOMER1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 125 amino acids of human homer homolog 1
Target Name	homer homolog 1
Other Names	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
Accession No.	Swissprot:Q86YM7Gene Accession:BC015502
Uniprot	Q86YM7
GeneID	9456;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only