

Recombinant Human HOMER2

Catalog No: #GP11631



Package Size: #GP11631-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HOMER2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 125 amino acids of human homer homolog 2
Target Name	homer homolog 2
Other Names	CPD; ACPD; VESL-2; HOMER-2
Accession No.	Swissprot:Q9NSB8Gene Accession:BC012109
Uniprot	Q9NSB8
GeneID	9455;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function. The encoded protein is a postsynaptic density scaffolding protein. Alternative splicing results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 14.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only