

Recombinant Human HSD17B7

Catalog No: #GP11638



Package Size: #GP11638-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HSD17B7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 180 amino acids of human hydroxysteroid (17-beta) dehydrogenase 7
Target Name	hydroxysteroid (17-beta) dehydrogenase 7
Other Names	PRAP; SDR37C1
Accession No.	Swissprot:P56937Gene Accession:BC007068
Uniprot	P56937
GeneID	51478;
Storage	-20~-80°C, pH 7.6 PBS

Background

HSD17B7 encodes an enzyme that functions both as a 17-beta-hydroxysteroid dehydrogenase in the biosynthesis of sex steroids and as a 3-ketosteroid reductase in the biosynthesis of cholesterol.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only