

Recombinant Human HSD3B1

Catalog No: #GP11645



Package Size: #GP11645-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HSD3B1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 280 amino acids of human hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 1
Target Name	hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 1
Other Names	I; HSD3B; HSDB3; HSDB3A; SDR11E1; 3BETAHSD
Accession No.	Swissprot:P14060Gene Accession:BC031999
Uniprot	P14060
GenID	3283;
Storage	-20~-80°C, pH 7.6 PBS

Background

3 beta-hydroxysteroid dehydrogenase (3-HSD), also known as HSD3B1 or HSDB3, is a bifunctional enzyme that plays a crucial role in the synthesis of all classes of hormonal steroids. Two human 3-HSD proteins, designated type I (3-HSD) and type II (3-HSD2), are expressed by different genes and function in different areas of the body.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only