

Recombinant Human TICAM2

Catalog No: #GP11653



Package Size: #GP11653-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TICAM2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 21-235 amino acids of human toll-like receptor adaptor molecule 2
Target Name	toll-like receptor adaptor molecule 2
Other Names	TIRP; TRAM; TIRAP3; MyD88-4; TICAM-2
Accession No.	Swissprot:Q86XR7Gene Accession:BC109266
Uniprot	Q86XR7
GenID	100302736;353376;
Storage	-20~-80°C, pH 7.6 PBS

Background

TIRP is a Toll/interleukin-1 receptor (IL1R; MIM 147810) (TIR) domain-containing adaptor protein involved in Toll receptor signaling

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only