

Recombinant Human ICAM4

Catalog No: #GP11654



Package Size: #GP11654-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ICAM4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-225 amino acids of human intercellular adhesion molecule 4
Target Name	intercellular adhesion molecule 4
Other Names	LW; CD242
Accession No.	Swissprot:Q14773Gene Accession:BC000046
Uniprot	Q14773
GenID	3386;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the Landsteiner-Wiener (LW) blood group antigen(s) that belongs to the immunoglobulin (Ig) superfamily, and that shares similarity with the intercellular adhesion molecule (ICAM) protein family. This ICAM protein contains 2 Ig-like C2-type domains and binds to the leukocyte adhesion LFA-1 protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only