

# Recombinant Human IRGM

Catalog No: #GP11677



Package Size: #GP11677-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human IRGM
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	immunity-related GTPase family, M
Other Names	IFI1; IRGM1; LRG47; LRG-47
Accession No.	Swissprot:A1A4Y4Gene Accession:BC128168
Uniprot	A1A4Y4
GeneID	345611;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the p47 immunity-related GTPase family. The encoded protein may play a role in the innate immune response by regulating autophagy formation in response to intracellular pathogens. Polymorphisms that affect the normal expression of this gene are associated with a susceptibility to Crohn's disease and tuberculosis.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only