

# Recombinant Human ITK

Catalog No: #GP11678



Package Size: #GP11678-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ITK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 170-340 amino acids of human IL2-inducible T-cell kinase
Target Name	IL2-inducible T-cell kinase
Other Names	EMT; LYK; LPFS1; PSCTK2
Accession No.	Swissprot:Q08881Gene Accession:BC109077
Uniprot	Q08881
GeneID	3702;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only