

Recombinant Human JAKMIP2

Catalog No: #GP11680



Package Size: #GP11680-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human JAKMIP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 215 amino acids of human janus kinase and microtubule interacting protein 2
Target Name	janus kinase and microtubule interacting protein 2
Other Names	NECC1; JAMIP2
Accession No.	Swissprot:Q96AA8Gene Accession:BC017354
Uniprot	Q96AA8
GeneID	9832;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only