

Recombinant Human LAMB3

Catalog No: #GP11706



Package Size: #GP11706-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LAMB3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human laminin, beta 3
Target Name	laminin, beta 3
Other Names	LAM5; LAMNB1; BM600-125KDA
Accession No.	Swissprot:Q13751Gene Accession:BC075838
Uniprot	Q13751
GeneID	3914;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product encoded by this gene is a laminin that belongs to a family of basement membrane proteins. This protein is a beta subunit laminin, which together with an alpha and a gamma subunit, forms laminin-5. Mutations in this gene cause epidermolysis bullosa junctional Herlitz type, and generalized atrophic benign epidermolysis bullosa, diseases that are characterized by blistering of the skin. Multiple alternatively spliced transcript variants that encode the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only