

Recombinant Human LCMT1

Catalog No: #GP11708



Package Size: #GP11708-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LCMT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 209 amino acids of human LCMT1
Target Name	leucine carboxyl methyltransferase 1
Other Names	LCMT; PPMT1; CGI-68
Accession No.	Swissprot:Q9UIC8Gene Accession:BC014217
Uniprot	Q9UIC8
GeneID	51451;
Storage	-20~-80°C, pH 7.6 PBS

Background

LCMT1 catalyzes the methylation of the carboxyl group of the C-terminal leucine residue (leu309) of the catalytic subunit of protein phosphatase-2A.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only