

Recombinant Human LDHC

Catalog No: #GP11709



Package Size: #GP11709-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LDHC
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 307 amino acids of human LDHC
Target Name	lactate dehydrogenase C
Other Names	CT32; LDH3; LDHX
Accession No.	Swissprot:P07864Gene Accession:BC064388
Uniprot	P07864
GeneID	3948;
Storage	-20~-80°C, pH 7.6 PBS

Background

Lactate dehydrogenase C catalyzes the conversion of L-lactate and NAD to pyruvate and NADH in the final step of anaerobic glycolysis. LDHC is testis-specific and belongs to the lactate dehydrogenase family. Two transcript variants have been detected which differ in the 5' untranslated region.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only