

# Recombinant Human LDHD

Catalog No: #GP11710



Package Size: #GP11710-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LDHD
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 62-265 amino acids of human LDHD
Target Name	lactate dehydrogenase D
Other Names	DLD
Accession No.	Swissprot:Q86WU2Gene Accession:BC047902
Uniprot	Q86WU2
GeneID	197257;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. The similar protein in yeast has both D-lactate and D-glycerate dehydrogenase activities. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only