

Recombinant Human LYN

Catalog No: #GP11727



Package Size: #GP11727-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 63-247 amino acids of human v-src-1 Yamaguchi sarcoma viral related oncogene homolog
Target Name	v-src-1 Yamaguchi sarcoma viral related oncogene homolog
Other Names	JTK8; p53Lyn; p56Lyn
Accession No.	Swissprot:P07948Gene Accession:BC126456
Uniprot	P07948
GeneID	4067;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only