

Recombinant Human MAGEB10

Catalog No: #GP11731



Package Size: #GP11731-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MAGEB10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 111-310 amino acids of human melanoma antigen family B, 10
Target Name	melanoma antigen family B, 10
Accession No.	Swissprot:Q96LZ2Gene Accession:BC101308
Uniprot	Q96LZ2
GeneID	139422;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only