

Recombinant Human BUB3

Catalog No: #GP11732



Package Size: #GP11732-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BUB3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 5-262 amino acids of human BUB3 mitotic checkpoint protein
Target Name	BUB3 mitotic checkpoint protein
Other Names	BUB3L; hBUB3
Accession No.	Swissprot:O43684Gene Accession:BC005138
Uniprot	O43684
GeneID	9184;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only