

Recombinant Human MAS1

Catalog No: #GP11737



Package Size: #GP11737-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MAS1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 246-325 amino acids of human MAS1 oncogene
Target Name	MAS1 oncogene
Other Names	MAS; MGRA
Accession No.	Swissprot:P04201Gene Accession:BC110454
Uniprot	P04201
GeneID	4142;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a class I seven-transmembrane G-protein-coupled receptor. The encoded protein is a receptor for angiotensin-(1-7) and preferentially couples to the Gq protein, activating the phospholipase C signaling pathway. The encoded protein may play a role in multiple processes including hypotension, smooth muscle relaxation and cardioprotection by mediating the effects of angiotensin-(1-7).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only