

Recombinant Human MATR3

Catalog No: #GP11738



Package Size: #GP11738-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MATR3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 250 amino acids of human matrin 3
Target Name	matrin 3
Other Names	MPD2; ALS21; VCPDM
Accession No.	Swissprot:P43243Gene Accession:BC015031
Uniprot	P43243
GeneID	9782;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a nuclear matrix protein, which is proposed to stabilize certain messenger RNA species. Mutations of this gene are associated with distal myopathy 2, which often includes vocal cord and pharyngeal weakness. Alternatively spliced transcript variants, including read-through transcripts composed of the upstream small nucleolar RNA host gene 4 (non-protein coding) and matrin 3 gene sequence, have been identified. Pseudogenes of this gene are located on chromosomes 1 and X.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only