

# Recombinant Human MBNL3

Catalog No: #GP11739



Package Size: #GP11739-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MBNL3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 100-250 amino acids of human muscleblind-like splicing regulator 3
Target Name	muscleblind-like splicing regulator 3
Other Names	CHCR; MBLX; MBXL; MBLX39
Accession No.	Swissprot:Q9NUK0Gene Accession:BC042090
Uniprot	Q9NUK0
GeneID	55796;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only