

## Recombinant Human OVAL

Catalog No: #GP11743



Package Size: #GP11743-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human OVAL
Brief Description	Recombinant Protein
Immunogen Description	OVAL protein
Target Name	ovalbumin (SERPINB14)
Other Names	SERPINB14
Accession No.	Swissprot:P01012Gene Accession:AAB59956
Uniprot	P01012
GeneID	396058;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor. It is presumed to be a storage protein.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only