

Recombinant Human NTMT1

Catalog No: #GP11749



Package Size: #GP11749-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NTMT1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	N-terminal Xaa-Pro-Lys N-methyltransferase 1
Other Names	NRMT; NTM1A; AD-003; HOMET1A; C9orf32; METTL11A
Accession No.	Swissprot:Q9BV86Gene Accession:BC001396
Uniprot	Q9BV86
GeneID	28989;
Storage	-20~-80°C, pH 7.6 PBS

Background

The METTL11A gene encodes an N-terminal methyltransferase for the RAN (MIM 601179) guanine nucleotide exchange factor regulator of chromosome condensation 1 (RCC1; MIM 179710). METTL11A enzyme alpha-N-methylates other protein targets such as SET.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only