

Recombinant Human MIER2

Catalog No: #GP11757



Package Size: #GP11757-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MIER2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 195-345 amino acids of human mesoderm induction early response 1, family member 2
Target Name	mesoderm induction early response 1, family member 2
Other Names	Mi-er2; KIAA1193
Accession No.	Swissprot:Q8N344Gene Accession:BC028203
Uniprot	Q8N344
GeneID	54531;
Storage	-20~-80°C, pH 7.6 PBS

Background

MIER2(mesoderm induction early response 1, family member 2), also known as Mi-er2 or KIAA1193, is a protein that in humans is encoded by the MIER2 gene. It is a transcriptional repressor.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only