

Recombinant Human MSR1

Catalog No: #GP11766



Package Size: #GP11766-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MSR1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 350-450 amino acids of human macrophage scavenger receptor 1
Target Name	macrophage scavenger receptor 1
Other Names	SRA; SR-A; CD204; phSR1; phSR2; SCARA1
Accession No.	Swissprot:P21757Gene Accession:BC063878
Uniprot	P21757
GenID	4481;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the class A macrophage scavenger receptors, which include three different types (1, 2, 3) generated by alternative splicing of this gene. These receptors or isoforms are macrophage-specific trimeric integral membrane glycoproteins and have been implicated in many macrophage-associated physiological and pathological processes including atherosclerosis, Alzheimer's disease, and host defense. The isoforms type 1 and type 2 are functional receptors and are able to mediate the endocytosis of modified low density lipoproteins (LDLs).?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only